



### Industrial Consultancy & Sponsored Research (IC&SR)

## QUADRATIC ELECTRO-OPTIC BASED DEFLECTION-FREE WIDE PATH-LENGTH MODULATION AND LATERAL SCANNING DEVICE FOR TIME DOMAIN OPTICAL COHERENCE TOMOGRAPHY

### IITM Technology Available for Licensing

#### Problem Statement

- Conventional systems involving FDOCT, Spectral Domain OCT (SDOCT) and Swept Source OCT (SSOCT), offers rapid axial scanning capabilities without mechanical components.
- However, it presents challenges such as **substantial data processing** requirements and the need for **high-resolution spectrometers** and line-scan cameras for **imaging deep tissues**, making the detection system bulky and complex

#### Technology Category/ Market

**Category**– Medical and Surgical/ Non Destructive Testing

**Applications** –Test Equipments, NDE, Biomedical systems, Sensors, Medical imaging

**Industry**- Biomedical Engineering, **Healthcare**

**Market** -The global medical imaging market size was valued at **USD 32.3 billion in 2022** and is expected to grow at a compound annual growth rate (CAGR) of **4.8% from 2023 to 2030**.

#### Key Features / Value Proposition

##### Technical Perspective:

- ❑ **A high speed axial and lateral scanning device** for time domain optical coherence(TCOCT) system that enables **non-mechanical noiseless imaging**
- ❑ **Interferometric imaging** provides **highly sensitive and can reveal fine details and subtle changes in the sample**, also allows for precise measurement of optical path differences and variations
- ❑ Involves **less data processing** and **does not have constraints on depth imaging** .
- ❑ **Polarized Imaging** uses **polarizers and quarter-wave plates** to manipulate the polarization state of the light beam that can reveal information about the sample's optical properties and anisotropy

##### User Perspective:

- ❑ Versatile technique and can be used for both biological samples that include **tissue analysis or cellular imaging as well as non-biological samples**.
- ❑ Offers **a range of imaging modalities**, making them valuable tools for scientific research, medical diagnostics, and materials characterization

#### Technology

- ❑ The present invention discloses electro-optic systems capable of performing axial, lateral, multi-dimensional, and even three-dimensional imaging
- ❑ The said electro-optic system for axial scanning of a sample comprising of:

**A light source** producing a light beam

A **detector** and an interferometer connected to the light source, where **Interferometer** includes a polarizer


A **beam splitter** that divides the light beam into two beams.

**Sample** defining a **sample arm**.

Electro-optic crystal maintained at a fixed temperature and linked to a **voltage source within the reference arm**

At least one **electro-optic crystal** is maintained at a predetermined temperature within the sample arm.

- ❑ Involves utilization of two sets of electro-optic crystals, where each set is maintained at two predetermined temperatures.
- ❑ The sample used is electro-optic crystal (6) is KTN crystal of the formula  $KTa_{1-x}Nb_xO_3$



**Fig.1** is the schematic representation of quadratic electro-optic based path-length modulation.

#### CONTACT US

**Dr. Dara Ajay, Head**

Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

**IITM TTO Website:**

<https://ipm.icsr.in/ipm/>

**Email:** [smipm-icsr@icsrpis.iitm.ac.in](mailto:smipm-icsr@icsrpis.iitm.ac.in)

[sm-marketing@imail.iitm.ac.in](mailto:sm-marketing@imail.iitm.ac.in)

**Phone:** +91-44-2257 9756/ 9719

Different types of imaging of imaging can be performed with the described electro-optic systems, include :

❖ **Axial Scanning Imaging:**

Allows for imaging along the axial (depth) direction of the sample and provides insights into the internal structure and composition of the sample.

❖ **Lateral Scanning Imaging:**

Enables imaging in two dimensions, typically the x and y axes useful for capturing cross-sectional or planar views of the sample's surface or structure.

❖ **Multi-Dimensional Imaging:**

Combines axial and lateral scanning, allowing for 3D imaging and provides a comprehensive view of the sample in terms of depth and lateral dimensions.

### Images




Fig.2 is the image of orientation and electro-optic phase tuning in KTN.




Fig.3(a) is the image of demineralised tooth sample




Fig.3(b) and 3(c) are images of demineralised tooth sample before mineralization and after 60hr demineralization B scans




Fig 4. Lateral screening and path length modulation in motionless two dimensional imaging




Fig.5 shows the schematic representation of electrooptically tuned lateral scanning




Fig.6 STD-OCT with quadratic electro-optic axial scanning

### Intellectual Property

- IITM IDF Ref. 1448
- IN380174-Granted

### TRL (Technology Readiness Level)

TRL-4, Technology Validated in the Lab

### Research Lab

Prof. Nilesh Jayantilal Vasa  
Department of Engineering Design  
Prof. Sarathi R  
Department of Electrical Engineering

### CONTACT US

Dr. Dara Ajay, Head  
Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:  
<https://ipm.icsr.in/ipm/>

Email: [smipm-icsr@icsrpis.iitm.ac.in](mailto:smipm-icsr@icsrpis.iitm.ac.in)

[sm-marketing@imail.iitm.ac.in](mailto:sm-marketing@imail.iitm.ac.in)

Phone: +91-44-2257 9756/ 9719