

### SYSTEM AND METHOD OF BEAMFORMING IN PULSE-ECHO ULTRASOUND IMAGING

#### IITM Technology Available for Licensing

#### Problem Statement

- Ultrasound imaging is crucial for diagnostics due to its safety, real-time imaging, and high resolution.
- Traditional Conventional Focused Transmit (CFT) has limited lateral resolution except at the focus.
- Synthetic Aperture (SA) techniques like **Synthetic Transmit Aperture (STA)** and **Multi-element Synthetic Transmit Aperture (MSTA)** improve lateral resolution but suffer from poor signal-to-noise ratio (SNR) and depth penetration.
- The problem is to explore Backprojection (BP) as an alternative beamforming method for SA-based ultrasound imaging to enhance SNR and depth while maintaining high resolution.

#### Intellectual Property

- IITM IDF Ref. **1595**
- **IN 380859 - Patent Granted**

#### Technology Category/ Market

**Category - Medical Imaging and Healthcare Applications** - Medical Imaging and Healthcare Industry - Medical Device Manufacturers


**Market** - The global ultrasound market has an estimated **CAGR of 6.3%** and a revenue size in the of \$8.5 billion in 2023.

#### TRL (Technology Readiness Level)

**TRL - 4: Technology validated in lab scale.**

#### Research Lab

**Prof. Arun K. Thittai,**  
Dept. of Applied Mechanics & Biomedical Engineering



**FIG.1. Shows an ultrasonic synthetic array aperture imaging system.**

#### Technology

The present invention discloses a pulse echo ultrasound imaging system and a method of beamforming in the pulse echo ultrasound imaging system.

For a Synthetic Transmit Aperture (STA) transmitting ultrasound diverging beams over a target, reflection data are received and stored in matrices for a plurality of transmit firings. (Fig. 1 & 2)

Image reconstruction is performed by weighting the backscattered echo signal and backprojecting on an ellipsoidal shell and averaging for every combination of transmit-receive pair.

The method implements an elliptical backprojection (EBP) formula that is derived from solving linear wave equation for STA.

For a Multi element Synthetic Transmit Aperture - Diverging beam (MSTA-DB) the eccentricity map is further modified with a virtual source placed behind the transducer. (Fig. 3)

#### CONTACT US

**Dr. Dara Ajay, Head**  
Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

**IITM TTO Website:**  
<https://ipm.icsr.in/ipm/>

**Email:** [smipm-icsr@icsrpis.iitm.ac.in](mailto:smipm-icsr@icsrpis.iitm.ac.in)  
[sm-marketing@imail.iitm.ac.in](mailto:sm-marketing@imail.iitm.ac.in)

**Phone:** +91-44-2257 9756/ 9719

### Key Features / Value Proposition

#### 1. Enhanced Image Quality:

The invention improves image quality by providing better Signal-to-Noise Ratio (SNR) and Contrast-to-Noise Ratio (CNR) compared to conventional methods, making it invaluable for accurate diagnostics.

#### 2. Superior Lateral Resolution:

The method offers superior lateral resolution, ensuring that detailed structures are clearly visualized, aiding in more precise medical assessments.

#### 3. Efficient Data Processing:

By utilizing matrices to store reflection data and implementing an elliptical backprojection (EBP) formula, the system optimizes data processing, enabling efficient and rapid image reconstruction.

#### 4. Adaptability:

The technique is adaptable for both Synthetic Transmit Aperture (STA) and Multi-element Synthetic Transmit Aperture - Diverging beam (MSTA-DB) systems, making it versatile for various ultrasound imaging applications.

#### 5. Competitive Advantage:

This innovation provides a competitive edge in the field of ultrasound imaging, potentially leading to improved patient care and better clinical outcomes.




FIG. 2. illustrates a schematic representation of the STA transmit, and a full aperture receive used in STA.




FIG. 3. Shows a schematic representation of the MSTA transmit, and full aperture receive used in MSTA.

### CONTACT US

Dr. Dara Ajay, Head

Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:

<https://ipm.icsr.in/ipm/>

Email: [smipm-icsr@icsrpis.iitm.ac.in](mailto:smipm-icsr@icsrpis.iitm.ac.in)

[sm-marketing@imail.iitm.ac.in](mailto:sm-marketing@imail.iitm.ac.in)

Phone: +91-44-2257 9756/ 9719