

### TRANSDERMAL COLLAGEN AND HEMOGLOBIN SENSOR

#### IITM Technology Available for Licensing

#### Problem Statement

- The diagnosis and assessment of skin disorders and the evaluation of cosmetic and anti-aging treatments **currently rely on invasive procedures like skin biopsies** or limited, expensive instruments with bulky instrumentation.
- Hence there is need for creating devices that enables **non-invasive, in-vivo detection** of disease biomarkers, such as **collagen and hemoglobin, in the subsurface layers of the skin.**

#### Technology Category/ Market

**Category – Medical and Surgical**

**Applications**–Biomedical systems, Cosmetics, Sensors, non-invasive techniques

**Industry – Healthcare, Biomedical devices**

**Market -Global Medical Devices Market size** was valued at USD 62.6 billion in 2021 and expected to grow from **USD 63.4 billion in 2022 to USD 134.56 billion by 2030**, at a **CAGR of 11.35% in the forecast period (2023-2030).**

#### Intellectual Property

- IITM IDF Ref. 1706
- IN201841024949

#### Research Lab

**Prof. SUJATHA N**

Dept. of Applied Mechanics

#### Key Features / Value Proposition

##### Technical Perspective

- ❑ Provides a transdermal fibre optic probe integrated reflectance collection device that is **non-invasive to differentiate normal and scleroderma skin**
- ❑ Spectral features of combined **reflectance and fluorescence spectrum** are compared with features of the training dataset for prediction of occurrence of the disease

##### User Perspective

- ❑ **Highly efficient and cost-effective** skin collagen and hemoglobin sensor device
- ❑ Capable of **assessing the skin condition** under various conditions of ailment during **changes in the subsurface collagen and hemoglobin content.**

#### Technology

The invention discloses an **apparatus for measuring/ determining the collagen and/or hemoglobin on a sample (Fig. 1); It includes:**

**Light Sources**-for exciting and illuminating sample

**Optical Probe** with a fiber optic probe and an integrated sphere

**Integrated Sphere**-collect reflected light (reflectance) from the sample.

**Sensor** -ensures that the probe is at the correct distance from the sample

**Spectrometer**-sense the collected fluorescence and reflected light

**Data Processing Unit** -processing data from the fluorescence and reflectance spectra.

**Display Unit**-displays in the form of fluorescence and reflectance spectra

- ❑ A **monochromatic light source** is used to **excite the sample**
- ❑ The Optical Probe consists of a **fiber optic probe and an integrated sphere.**
- ❑ The fiber probe is connected to the monochromatic light source and has two main sections- **Center Section** having fiber cable that connects to the excitation source for the sample and **Outer Section that** contains a bundle of fiber cables arranged in a circular pattern at the periphery of the center section.




Fig. 1

#### CONTACT US

**Dr. Dara Ajay, Head**

Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

**IITM TTO Website:**

<https://ipm.icsr.in/ipm/>

**Email:** [smipm-icsr@icsrpis.iitm.ac.in](mailto:smipm-icsr@icsrpis.iitm.ac.in)

[sm-marketing@imail.iitm.ac.in](mailto:sm-marketing@imail.iitm.ac.in)

**Phone:** +91-44-2257 9756/ 9719

## Industrial Consultancy & Sponsored Research (IC&SR)

Further, discloses a method of operation for fluorescence and reflectance measurement that includes:

### Fluorescence Measurement

- Positioning the sensor
- Excitation the sample
- Collect the emitted fluorescence light
- Direct the fluorescence light to the spectrometer
- Processing data and Display


### Reflectance Measurement

- Illuminate the sample using the light source
- Collect the reflected light from the sample using the integrated sphere
- Directing the collected reflected light into the spectrometer.
- Process the data and Display the reflectance spectrum


### Analysis

- Compare and analyze the spectral features of the reflectance and fluorescence spectra
- Utilizing a pre-loaded algorithm for data processing within the processing unit.
- Display the results


### Images





**Fig.2:** illustrates acquisition of fluorescence spectrum from (a) dorsal palm (b) forearm and (c) above the elbow



**Fig. 3** illustrates acquisition of diffuse reflectance spectrum from (a) dorsal palm (b) forearm and (c) above the elbow



**Fig. 4:** illustrates Stained microscope images obtained from abnormal Fig 4 (a); and normal biopsy site Fig. 4(b)



**Fig. 5:** illustrates (a) fluorescence mode and (b) DRS mode

### TRL (Technology Readiness Level)

**TRL-5, Technology Validated in Relevant Environment**

### CONTACT US

**Dr. Dara Ajay, Head**  
Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

**IITM TTO Website:**  
<https://ipm.icsr.in/ipm/>

**Email:** [smipm-icsr@icsrpis.iitm.ac.in](mailto:smipm-icsr@icsrpis.iitm.ac.in)

[sm-marketing@imail.iitm.ac.in](mailto:sm-marketing@imail.iitm.ac.in)

**Phone:** +91-44-2257 9756/ 9719