

A PROCESS FOR FABRICATION OF PURE/MACROPOROUS APATITIC (CDHA) BONE CEMENT FOR NON-LOAD BEARING ORTHOPEDIC APPLICATIONS

IITM Technology Available for Licensing

Problem Statement

- Addressing bone tissue defects and loss is a significant concern in orthopedics, requiring advanced regenerative procedures.
- Highlighting the promising features of CPCs, such as higher biocompatibility, bioactivity, and bioresorbability, for bone grafting applications.
- Exploring the potential of utilizing eggshell waste as a rich source of calcium for deriving calcium phosphate materials with trace elements mimicking human bone composition.
- Further, Recognizing limitations in current bone cement technologies and proposing a novel process for fabricating calcium deficient hydroxyapatite (CDHA) bone cements with improved injectability, macroporosity, and resorbability properties for effective bone grafting applications.

Technology Category/ Market

Biomaterials and Biomedical Engineering

Applications- Bone Grafting, Implant Fixations, Void Filling, Orthopedic Coating

Industry - Orthopedic Devices and Implants.

Market- Global orthopedic biomaterials market size was estimated at USD 19.2 billion in 2022 and is expected to grow at a **CAGR of 7.8%** from 2023 to 2030.

Research Lab

Prof. Sampath Kumar T S

Dept. of Metallurgical and Materials Engineering

Prof. Mukesh Doble

Dept. of Biotechnology




Fig. 1. Represents (a) Syringe containing cement paste, (b) Extruded injectable bone cement paste and (c) Cylindrical and round shaped molds of synthetic bone cement.

Intellectual Property

- IITM IDF Ref. 2276
- IN 410938 - Patent Granted

TRL (Technology Readiness Level)

TRL - 4: Technology validated in lab scale.

Technology

- The present invention provides a **process for fabrication of macroporous apatitic (CDHA) bone cement** comprising steps of:
 - (i) Synthesizing Nanocrystalline Hydroxyapatite (HA) using calcium nitrate ($\text{Ca}_3(\text{NO}_4)_2 \cdot 4\text{H}_2\text{O}$) and di-ammonium hydrogen phosphate ($(\text{NH}_4)_2\text{HPO}_4$) which are mixed and subjected to irradiation in a microwave oven of 800 W for 30 minutes. Thereafter, washing the precipitate with distilled water followed by oven dried at 100 °C overnight and powdered finely.

CONTACT US

Dr. Dara Ajay, Head

Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:
<https://ipm.icsr.in/ipm/>

Email: smipm-icsr@icsrpis.iitm.ac.in

sm-marketing@imail.iitm.ac.in

Phone: +91-44-2257 9756/ 9719



Technology (Contd.)

- (ii) Synthesizing pure eggshell derived tricalcium phosphate (ESD β -TCP) by heating a mixture of CaO extracted from eggshell and anhydrous dicalcium phosphate in 1:2 molar ratio at 1000 °C for 12 hours.
- (iii) Preparing a workable dough of the cements by mixing solid phase (β -TCP and HA) and liquid phase with a porogen agent, wherein cements (after homogenous mixing) were allowed to set for 1 hour and then, incubated in the PBS (pH ~7.4) solution at physiological conditions for certain time intervals; thereafter, the cements were processed for its further characterization and biocompatibilities investigations.

Key Features / Value Proposition

•The process introduces a microwave-accelerated wet chemical synthesis for nano-crystalline hydroxyapatite, ensuring efficient and controlled production.

**Innovative
Synthesis
Methodology**

•Utilizing eggshell waste for synthesizing pure β -tricalcium phosphate (ESD β -TCP) adds sustainability and cost-effectiveness to the manufacturing process.

**Eggshell-Derived
 β -TCP Synthesis**

•The cement preparation involves a precisely formulated combination of synthetic and eggshell-derived components, providing versatility and adaptability for varied orthopedic applications.

**Tailored Cement
Formulation**

•Incorporating disodium hydrogen phosphate, gelatin, and chitosan in the liquid phase enhances the mechanical and biological properties of the resulting bone cement.

**Polymeric Liquid
Phase
Enhancement**

•The inclusion of mannitol or polysorbate as a porogen agent introduces a controlled porosity element, contributing to improved macroporosity in the final product.

**Innovative
Porogen Agent
Integration**

•The process ensures the formation of a pure apatitic (CDHA) phase in bone cement with defined dimensions, allowing for efficient bioresorption and enhanced biocompatibility.

**Optimized CDHA
Phase
Formation**

CONTACT US

Dr. Dara Ajay, Head

Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:

<https://ipm.icsr.in/ipm/>

Email: smipm-icsr@icsrpis.iitm.ac.in

sm-marketing@iimail.iitm.ac.in

Phone: +91-44-2257 9756/ 9719