

### MULTIFUNCTIONAL HYDROGEL FOR WOUND HEALING

#### IITM Technology Available for Licensing

#### Problem Statement

- Current wound care lacks optimal solutions, with traditional dressings often ineffective for proper healing, necessitating advancements in wound care technology.
- Injectable hydrogels show promise for wound dressings but suffer from limitations such as fragility and complexity in composition, hindering widespread adoption and efficacy.
- There is a need for a **pH-responsive self-healing, injectable, multifunctional hydrogel** that addresses the complexities of wound healing.

#### Intellectual Property

- IITM IDF Ref. 2780
- IN 202441004292

#### Technology Category/ Market

##### Category - Advanced Wound Care

**Applications** - Chronic Wound Management, Acute Wound Healing, Dermatology and Cosmetic Applications.

**Industry** - Healthcare and Medical Devices, Pharmaceuticals and Biotechnology


**Market**- Advance wound care market is expected to be valued at US\$ 14,667.7 million by 2034 with a growth at a **CAGR of 2.4%**.

#### TRL (Technology Readiness Level)

**TRL-4, Technology validated in relevant environment.**

#### Research Lab

**Prof. Vignesh Muthuvijayan**  
Dept. of Biotechnology



**FIG. 1** illustrates self-healing potential of hydrogel (A) Photograph showing the self-healing efficiency of the hydrogels. (B)  $G'$  and  $G''$  of the CS-o-CMC hydrogels when alternate step strain switched from a small strain ( $\gamma = 10\%$ ) to a large strain ( $\gamma = 100\%$ ) at a fixed angular frequency (10 rad/s).

#### Technology

- The present invention pertains to a **pH-responsive, self-healing and injectable hydrogel for wound healing**, comprising a combination of polysaccharides and a biomolecule, specifically:
  - **Chitosan**, facilitating cell proliferation support;
  - **Oxidized carboxymethyl cellulose**, enhancing rheological behavior; and
  - **Lauric acid**, imparting immunomodulatory properties.
- The hydrogel demonstrates self-healing capabilities, achieved through the formation of Schiff base and electrostatic interaction.
- The Schiff base bonds resulting from a chemical reaction between the amine of chitosan and the aldehyde of oxidized carboxymethyl cellulose.
- The electrostatic interaction forms between the free amine group of chitosan and  $\beta$ -glycerol phosphate.

#### CONTACT US

**Dr. Dara Ajay, Head**  
Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

**IITM TTO Website:**  
<https://ipm.icsr.in/ipm/>

**Email:** [smipm-icsr@icsrps.iitm.ac.in](mailto:smipm-icsr@icsrps.iitm.ac.in)  
[sm-marketing@iitm.ac.in](mailto:sm-marketing@iitm.ac.in)

**Phone:** +91-44-2257 9756/ 9719

### Method

1

- Uniformly mixing chitosan, oxidized carboxymethyl cellulose, and lauric acid using luer-lock syringe;

2

- Allowing the formation of Schiff base bonds through a reaction between the amine of chitosan and the aldehyde of oxidized carboxymethyl cellulose;

3

- Allowing the formation of electrostatic interaction at 37 °C using a cross-linker, specifically 40-50%  $\beta$ -glycerol phosphate to enhance stability and gelation time; and

4

- Utilizing a solvent, including MilliQ water or PBS or a combination of both, for dissolving the polymers.

### Key Features / Value Proposition

#### 1. Versatile Cross-Linking:

Double cross-linking such as Schiff base and  $\beta$ -glycerol phosphate offers versatility and stability of hydrogel in different wound environment. Importantly, pH of the formulation is intrinsically raised to 7.0-7.4.

#### 2. pH Responsive self-healing:

Alkaline pH-responsive self-healing provides efficient scaffolding system for the migration of cells and ensures the stability of hydrogel without breaking during wear and tear condition at the wound site.

#### 3. Controlled Drug Release:

Controlled release of immunomodulatory molecule in pH 8.5, aids in inflammation control and immune response modulation

#### 4. Enhanced Healing:

Injectable, self-healing hydrogel formulation accelerates wound healing by reducing inflammation, which promotes granulation tissue maturation, and re-epithelialization.




FIG. 1 (C) Diffusion of dynamic covalent imine bond in different pH medium at different time interval.




FIG. 2 illustrates in vivo wound healing potential of hydrogel. Photographs of wound contraction upon treatment with hydrogels at different time interval.

### CONTACT US

Dr. Dara Ajay, Head

Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:

<https://ipm.icsr.in/ipm/>

Email: [smipm-icsr@icsrpis.iitm.ac.in](mailto:smipm-icsr@icsrpis.iitm.ac.in)

[sm-marketing@imail.iitm.ac.in](mailto:sm-marketing@imail.iitm.ac.in)

Phone: +91-44-2257 9756/ 9719