

### METHOD FOR BLOOD PLASMA SEPARATION BASED ON ACOUSTOCAPILLARY AND ASYMMETRIC CAPILLARY FLOW IITM Technology Available for Licensing

#### Problem Statement

- Generally, there is major interest to generate high quality plasma from whole blood in many biomedical analysis & clinical diagnostic methods.
- In conventional method of centrifuge-based separation of plasma, there is commonly needs milliliters of blood & a **labor-intensive** handling process in addition to the **bulky apparatus**, which makes the process **difficult & unsustainable** for regular tests & thus sample pre-treatment step becomes a **bottleneck** of the assay process.
- Hence, there is a need to address above issues in efficient manner.

#### Technology Category/ Market

**Technology:** Microfluidic Device;

**Industry:** Pharmaceutical, Medical device;

**Application:** Medical Test Equipment;

**Market:** The global **microfluidic device** market is projected to reach **\$158.1B** by **2031**, at a **CAGR** of **22.4%** during (2023-2031)

#### Technology

- Present invention describes a **microfluidic device for separation of plasma from blood cells**.
- Said Patent further discloses a **method for blood plasma separation** using microfluidic device based on **acoustocapillary and asymmetric capillary flow**.
- Said **microfluidic device** comprises a **polydimethylsiloxane (PDMS) microchannel layer** bonded with a **PDMS coated glass** thereby making the glass walls hydrophilic slide and the bottom wall of the microchannel is hydrophobic.
- The **acoustic radiation force** on the blood cells act towards the center of the channel or node makes blood cells to concentrate at the center of the channel & cell-free plasma at the-

walls or anti-node and the acoustic waves propagate in the width wise direction which makes the cells to move in the y direction.

- The **method for blood plasma separation** using microfluidic device is illustrated in below smart chart & figures:





Fig 1

#### Intellectual Property

IITM IDF Ref. 1385; Patent No:410407;

TRL (Technology Readiness Level)

TRL-3, Proof of concept tested in Lab;

Research Lab

**Prof. Ashis Kumar Sen,**  
Dept. of Mechanical Engineering

#### CONTACT US

**Dr. Dara Ajay, Head TTO**  
Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

**IITM TTO Website:**  
<https://ipm.icsr.in/ipm/>


Email: [smipm-icsr@icsrpis.iitm.ac.in](mailto:smipm-icsr@icsrpis.iitm.ac.in)  
[sm-marketing@imail.iitm.ac.in](mailto:sm-marketing@imail.iitm.ac.in)  
Phone: +91-44-2257 9756/ 9719

### Key Features / Value Proposition


#### ❖ *Technical Perspective & Industrial Perspective:*

- Present invention facilitates **simple diagnostic process cum handy device** including features like **improved quality, reproducibility & reliability** of the assay results. (Refer Fig.1 illustrates acousto-capillary based blood plasma separation process.)
- Provides a capillary driven blood plasma separation device designed both in **horizontal configuration & vertical configuration**.
- **Enabled realization of a lab on chip diagnostic platform.**
- The cell-free plasma at the walls can be separated from the blood cells using a **trifurcated channel configuration &** separated plasma can be collected at the **outlets** for further analysis or detection.
- A **suction chamber** can also be used to **control fluid flow** along with the **capillary flow**.
- **Cost-effective & time-consuming process.**

### Image



**Fig. 2 : Illustrates asymmetric capillary flow based blood plasma separation in horizontal configuration.**



**FIGURE 3**

**Fig. 3: Illustrates asymmetric capillary flow-based blood plasma separation in vertical configuration.**

### CONTACT US

Dr. Dara Ajay, Head TTO  
Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

**IITM TTO Website:**  
<https://ipm.icsr.in/ipm/>

Email: [smipm-icsr@icsrpis.iitm.ac.in](mailto:smipm-icsr@icsrpis.iitm.ac.in)  
[sm-marketing@imail.iitm.ac.in](mailto:sm-marketing@imail.iitm.ac.in)  
Phone: +91-44-2257 9756/ 9719