



IMPROVED PERFUSION CUM COMPRESSION BIOREACTOR SYSTEM FOR TISSUE ENGINEERING APPLICATIONS

IITM Technology Available for Licensing

PROBLEM STATEMENT

- **Bio-reactors** are systems where conditions are **closely controlled** to permit and induce a certain behavior in **living cells or tissues**.
- The system is provided with **controlled and steady flow** of growth media and factors necessary for cell growth where the **pH, temperature, pressure, nutrient supply and waste removal** is maintained at optimum conditions.
- Some of the **types of bioreactors** available are spinner flasks, rotating vessels, hollow fiber, perfusion bioreactors .
- The present invention is an **improved design of a perfusion cum compression bio-reactor** for creating tissue like **cartilage tissues**.

TECHNOLOGY CATEGORY MARKET

Category: Biotechnology & Genetic Engineering

Industry: Biomedical devices Manufacturing, Cell-Based Therapy Manufacturing, Vaccine Manufacturing

Application: Tissue engineering, Biomedical engineering applications and basic science studies.

Market: The global market size of Biotechnology & Genetic Engineering was valued at **USD 1.36 Billion in 2023** and projected to grow from **USD 1.68 Billion in 2024 to USD 7.73 Billion by 2032**, exhibiting a **(CAGR) of 20.94%** during the forecast period (2024 - 2032).

INTELLECTUAL PROPERTY

IITM IDF Ref. 1556; Patent No: IN 532512;

TRL (Technology Readiness Level)

TRL-4, Experimentally validated in Lab;

Research Lab

Dr. Merlin Rajesh Lal L P ,
Dept. of Biotechnology, IIT Madras.

TECHNOLOGY



Figure 1 shows a **Assembled PCC Bioreactor chamber**

❖ **Perfusion Cum Compression (PCC)** bioreactor chamber was designed and fabricated with the following parts.

- 1) **Compression chamber,**
- 2) **Inflatable membrane O-ring,**
- 3) **Compression port, and**
- 4) **Lock plate.**

1. Compression chamber




Figure 2 shows a compression chamber

- The **acrylic** compression chamber was **slotted** with **12 number** of grew in **two steps**.
- First slot have a larger diameter and the lower with smaller diameter which hold **cell seeded scaffolds/tissue explants**.
- The **first step** has a cylindrical slot with diameter of **21mm** and the **second step** has a diameter of **610mm** **Perpendicular** to the lower grew the orifice of medium perfusion was placed.
- The both ends of the orifice (towards and from the cell seeded scaffold/tissue explants side) **the curvature is increased** to make the **flow of medium uniform** and outer end is enabled for threaded luer connectors with **silicon O-rings**.
- Silicon rubber tubing connected to the luer of the chamber pass through a **peristaltic pump** and then to the growth media reservoir.

CONTACT US

Dr. Dara Ajay, Head TTO
Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:
<https://ipm.icsr.in/ipm/>

Email: smipm-icsr@icsrpis.iitm.ac.in

sm-marketing@imail.iitm.ac.in

Phone: +91-44-2257 9756/ 9719

2. Inflatable membrane O-ring



Figure 3 shows a Inflatable membrane O-ring

- ❑ The unique inflatable membrane **O-ring** is fabricated using **silicon-rubber**.
- ❑ The **one end is inflatable** with **compressed air** and the other provisions to act as O-ring which gets suited to **the larger surface diameter** slot in the compression chamber.
- ❑ The dimensions of the slots **vary with the cell seeded** scaffolds/tissue explants.

3. Compression port




Figure 4 shows a Compression port

- 103- Compression port
- 17- Orifice of Compression port
- 18-Provision for connecting (barb) tubing

- The compression port (20mm Dia) was designed like an **inverted half head candle holder** made of **polypropylene or teflon or acrylic**.
- The broad end of the compression port gets suited **above** the Inflatable membrane O-ring.
- On the other end of the compression port, provisions are made for **connecting (barb) tubing** that supply compressed air.
- O-ring through the centre an **orifice of 3 mm** diameter was provisioned for supplying of compressed air.

4. Lock plate





Figure 5 shows a Lock Plate

- 104-Lock plate
- 19-Access for compression port
- 20-Fasteners hole on plate

- ❑ **Lock plate (2mm thick)** is made of **stainless steel** to lock the bioreactor chamber assembly intact.
- ❑ **Two plates, one at the bottom** of the compression chamber and the **other was suited above**.
- ❑ Using stainless steel fasteners to lock plates were locked keep the assembly intact.

Key Features / Value Proposition

- ✓ Can be used for **biomechanical forces** like **dynamic or static compression loading, hydrostatic pressure**.
- ✓ Growth medium of perfusion may be **separately or in combination** on **cell seeded scaffolds/tissue explants** for tissue engineering applications.
- ✓ The **compression chamber of size ratio 30x60x180** is **slotted** with **12 cylindrical slots**.
- ✓ **Minimum loading time** of the cell seeded scaffolds **less than** that **5-10 minutes**.

- ✓ **Compact** and consists of members which have dual functions.
- ✓ **vivo microenvironment** by **controlling the amount and type of biomechanical forces**.
- ✓ **Improved** bio-reactor for generating tissues like **neo-cartilaginous tissues, bone tissue, liver tissues** etc.
- ✓ **Easy to assemble, disassemble**.

CONTACT US

Dr. Dara Ajay, Head TTO
Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:
<https://ipm.icsr.in/ipm/>

Email: smipm-icsr@icsrpis.iitm.ac.in

sm-marketing@imail.iitm.ac.in

Phone: +91-44-2257 9756/ 9719