



IIT MADRAS

Indian Institute of Technology Madras

Technology Transfer Office TTO - IPM Cell



Industrial Consultancy & Sponsored Research (IC&SR)

SYNTHESIS OF TOXIC-FREE STAR SHAPED GOLD NANOSTRUCTURES USING MICROFLUIDICS AND ITS USAGE IN INTRACELLULAR DELIVERY

IITM Technology Available for Licensing

PROBLEM STATEMENT

- **Anisotropic gold nanoparticles**, with near-infrared tunability, are useful in cellular transfection and biomedical applications, **targeting cancer cells** and causing thermomechanical damage through laser radiation.
- **Gold nanoparticles** are created through **chemical reduction using reducing agents** like trisodium citrate, sodium borohydride, or acetic acid, **but controlling kinetics and producing nanostars** is challenging.
- **Microfluidics** technology offers **high yields** and low dispersion, **but** it requires large volumes of reagents and time.
- **A method to synthesize toxic-free star-shaped gold nanostructures is needed**, requiring **fewer reagents, reproducibility, cost-effectiveness, and time efficiency.**

TECHNOLOGY CATEGORY MARKET

Technology: Toxic-free star shaped gold nanostructures

Category: Micro & Nano Technologies

Industry: Bio-Micro/Nano, Biomedical Engineering

Application: Intracellular delivery using photoporation

Market: The global market size of **nanomaterials in personalized medicine** is expected to **increase from \$346.7 billion in 2023 to \$592.0 billion** by the end of **2028**, with a compound annual growth rate (CAGR) of **11.3%** during the forecast period of **2023-2028**.

INTELLECTUAL PROPERTY

IITM IDF Ref. 2271, Patent No: IN 526929

TRL (Technology Readiness Level)

TRL-3, Experimental proof of concept ;

Research Lab

Dr. Tuhin Subhra Santra,
Dept. of Engineering Design.

TECHNOLOGY

A **method to synthesize toxic-free star shaped gold nanostructures** comprising steps of

1 **Preparing seed nanoparticles** by adding a reducing agent into a gold precursor;

2 **Diluting the seed solution** prepared in step (1);

3 **Adding the diluted seed solution** prepared in step (2) to a first inlet of droplet microfluidic device;

4 **Maintaining a low pH** throughout the synthesis by feeding a strong acid through a second inlet.

5 **Adding gold precursor** through a third inlet;

6 **Simultaneously adding ascorbic acid and silver nitrate in the next junction** through a fourth and fifth inlet, respectively;

7 **Adding silicone oil to the mixed reagents solution** and forming aqueous droplets in oil;

8 **To halt the growth process**, collecting the nanostar solution in PEG solution kept at ice-cold temperature; and

9 **Removing excess oil and PEG** by multiple centrifugations

CONTACT US

Dr. Dara Ajay, Head TTO
Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:
<https://ipm.icsr.in/ipm/>

Email: smipm-icsr@icsrpis.iitm.ac.in

sm-marketing@imail.iitm.ac.in

Phone: +91-44-2257 9756/ 9719

Microfluidic Device




Figure 1 illustrates synthesis of toxic-free gold nanostars (Au NS) in PDMS-based droplet microfluidic device




Figure 2 illustrates characterization of synthesized Au NS: (a) and (b) are transmission electron microscopy (TEM) images and enlarged portion of a thorn of single NS




Figure 3 illustrates characterization of synthesized Au NS (c) UV-Vis absorption spectrum and (d) DLS size distribution (lognormal) of NS.

Key Features / Value Proposition

❖ Application

- ✓ Effectively useful for **intracellular delivery**

❖ Toxic Content

- ✓ Toxic-free

❖ Shape

- ✓ **Star shaped** gold nanostructures and **crystalline nature**

❖ Deliver mode

- ✓ Deliver biomolecules to the cells by **light activated technique** mediated by star shaped gold nanostructures

❖ Synthesis

- ✓ **Droplet based** microfluidic device can be used to synthesize star shaped gold nanoparticles in a single step **without the use of any surfactant**

❖ Gold precursor

- ✓ **Tetrachloroauric acid**

❖ Strong acid is hydrochloric acid

❖ Thiolated polymer

- ✓ Poly(ethylene glycol) thiol

❖ Reducing agent

- ✓ **Trisodium citrate** for **Au seed formation**
- ✓ **Acetic acid** for **Nanostar formation**

❖ Flow rates

- ✓ **150 to 900 $\mu\text{L}/\text{h}$** to achieve the **oil to aqueous phase** flow rate ratio as 1 to 6.
- ✓ **At any point** of synthesis, parameters **can be tuned** by varying the flow rate of reagents and concentration ratio.

❖ Time efficiency

- ✓ Au NS synthesis in microfluidic device takes **only few seconds & cost-effectiveness.**

CONTACT US

Dr. Dara Ajay, Head TTO

Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:

<https://ipm.icsr.in/ipm/>

Email: smipm-icsr@icsrpis.iitm.ac.in

sm-marketing@imail.iitm.ac.in

Phone: +91-44-2257 9756/ 9719