



Industrial Consultancy & Sponsored Research (IC&SR)

2D SUBSTRATE PROTEIN MICROPATTERNING PROCESS FOR CELL ALIGNMENT

IITM Technology Available for Licensing

PROBLEM STATEMENT

- Cells interact with the **Extracellular Matrix (ECM)** through **receptors like Integrin**, sensing mechanical and chemical stimuli.
- **Organs have unique structures and alignments**, making **2D and 3D patterned microenvironments crucial for cell growth** in tissue-engineered organ substitutes.
- **Tissue engineering techniques** create customized microenvironments for **cell alignment**, varying from organ to organ.
- **Techniques** include electromagnetic fields, electrophoretic force, nanofiber scaffolds, 3D gels, and microstructure substrates.
- There is **no techniques using chemical stimulus** to achieve cell affinity and align cells on a substrate.
- The **need for an improved 2D substrate protein micropatterning process** exists to **align cells in a specific pattern** by providing chemical stimulus, achieving cell affinity to the specific stimulus.

TECHNOLOGY CATEGORY MARKET

Technology: Process for fabricating a SU 3D structure on a silicon substrate

Category: Precise single cell analysis

Industry: Biomedical industry

Application: Tissue engineering ,nanotechnology

Market: The global market size was estimated to be worth **\$3.5 billion in 2023** and is **poised to reach \$7.1 billion by 2028**, growing at a **CAGR of 15.3% from 2023 to 2028**

INTELLECTUAL PROPERTY

IITM IDF Ref. 1966

Patent No: IN 486480

TRL (Technology Readiness Level)

TRL- 4, Experimentally validated in Lab;

Research Lab

Prof. Tuhin Subhra Santra,
Dept. of Engineering Design

TECHNOLOGY

2D Substrate Protein Micropatterning Process

1

•Forming mold using SU8 3D structure on silicon substrate.

2

•Casting polydimethylsiloxane (PDMS) stamps using SU8 3D structure.

3

•Cleaning silicon substrate with Piranha solution at 80°C for 10 minutes and washing in deionized water.

4

•Dried with nitrogen blow, kept at 120 °C for 5 minutes for dehydration in order to form the SU8 3D structure

5


•Casting PDMS pillar array using mold.

6

•Aligning cells in a specific pattern using 2D substrate.

7

•Achieving cell affinity to specific chemical stimulus



The figure above illustrates a graphical representation of a final 3D SU8 mold on silicon substrate showing an array of through-holes in SU8 observed under optical microscope

CONTACT US

Dr. Dara Ajay, Head TTO

Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:


<https://ipm.icsr.in/ipm/>

Email: smipm-icsr@icsrpis.iitm.ac.in

sm-marketing@imail.iitm.ac.in


Phone: +91-44-2257 9756/ 9719

Fluorescence image of the pattern printed by the PDMS stamp on a glass substrate



- **SU8 3005** (Micro chemical) selected for 3D mold creation.
- Spin-coated on silicon sample at **1000 rpm for 60 seconds**.
- Soft-baked at **95 °C for 3 minutes**.
- Patterned with a mask using **UV lithography machine**.
- Hard-baked at **65 °C for 1 minute**, followed by 95 °C for 2 minutes.
- Cooled to room temperature and developed in **SU8 developer for 3 minutes**.
- Substrate washed with Isopropanol and **dried with nitrogen blow**.
- Optical microscope observed final sample with through-holes.

DMS stamping procedure for cell patterning



Mixing PDMS with a curing agent and

Treating SU8 substrate with trichlorosilane.

Transferring the mixture to the substrate.

Testing the final stamp with fluorescein ink.

Activating the stamp with protein solution.

Printing the protein on a tissue culture plate.

Seeding Siha cells on the patterned plate.

Staining cells with Calcein AM and observing under a fluorescence microscope.

Key Features / Value Proposition

□ Structure

- A mold is formed using a SU8 3D structure on a **silicon substrate** for casting polydimethylsiloxane (**PDMS stamps**).
- The **SU8 3005** (Micro chemical) is selected for making the 3D mold.
- The **SU8 3D** structure is **spin-coated, soft-baked, patterned, hard-baked**, and developed in a SU8 developer.

□ Cell Alignment

- Innovative Cell Alignment Process in **Extra Cellular Matrix (ECM)** Structure
- Improved 2D substrate protein micropatterning process for cell alignment.

□ Cell affinity

- Provides **chemical stimulus** to cells, enabling them to align in a specific pattern and achieve cell affinity.

□ Toxicity

- The process is **nontoxic**, easy to implement, and **can pattern heterogeneous cell lines**.

□ Performance

- Efficient analysis of **genetic disorder for the patient at a single cell level**
- Precise **genome and transcriptome analysis** at the single cell level, provide **functional consequence of mutation** and copy number variation of cells

CONTACT US

Dr. Dara Ajay, Head TTO
Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:
<https://ipm.icsr.in/ipm/>

Email: smipm-icsr@icsrpis.iitm.ac.in

sm-marketing@imail.iitm.ac.in

Phone: +91-44-2257 9756/ 9719