

# TTO - IPM Cell



# Industrial Consultancy & Sponsored Research (IC&SR)

# SYNTHESIS OF PROTEIN PROTECTED LUMINESCENT METAL CLUSTERS AND RETAINING THE BIOACTIVITY OF THE PROTEIN

IITM Technology Available for Licensing

## **PROBLEMSTATEMENT**

- Protein-protected noble metal nanoclusters (NCs) are a new class of biohybrid materials due to their intense luminescence, high quantum yield, stability, and biocompatibility.
- NCs have potential applications in sensing, biology, and multimodal imaging.
- > Protein-directed cluster synthesis has been explored for therapeutic applications, targeted biolabeling, and protein-protected clusters.
- Large proteins with high molecular weight, small proteins with low molecular weight, and gold are commonly used.
- Bioactivity loss due to irreversible conformational changes limits their widespread use in biomedical research.
- > A new class of atomically precise silver nanoclusters co-protected by hydride and phosphines as ligands has been reported

### TECHNOLOGYCATEGORY MARKET

**Technology:** Internalization of silver cluster by proteins and luminescent retaining bioactivity Category: Chemistry & Chemical analysis

Industry: Chemical/Pharmaceutical

**Application:** Biomedical

Market: The global market size was valued at USD 1.55 trillion in 2023 and is projected to grow at a compound annual growth rate (CAGR) of 13.96% from 2024 to 2030.

## **INIELLECTUAL PROPERTY**

IITM IDF Ref. 1779 **Patent No: IN 495343** 

TRL (Technology Readiness Level)

**TRL- 4**, Experimentally validated in Lab;

# Research Lab

Prof. Pradeep T. Dept. of Chemistry

### TECHNOLOGY

A method of synthesizing protein protected luminescent metal clusters, the said method comprises:

a preformed ligand protected atomically precise metal cluster as a precursor;


a protein, mixed with the metal cluster precursor

at pH 7.1 under continuous magnetic stirring for 12h at room temperature;

wherein the reaction mixture gradually change the color during the course of reaction,accompanied by the evolution of red luminescence of the solution upon continuous stirring,

which leads to the formation of red luminescent metal cluster indicating that the cluster core is incorporated in a single protein moiety without altering the secondary structure of the protein and its bioactivity

#### Schematic of the formation of Ag18@BSA



# **CONTACT US**

Dr. Dara Ajay, Head TTO Technology Transfer Office, IPM Cell- IC&SR, IIT Madras **IITM TTO Website:** 

https://ipm.icsr.in/ipm/

Email: headtto-icsr@icsrpis.iitm.ac.in

tto-mktg@icsrpis.iitm.ac.in


Phone: +91-44-2257 9756/ 9719



# TTO - IPM Cell



# Industrial Consultancy & Sponsored Research (IC&SR)



# Key Features / Value Proposition

## **Protein Bioactivity and Cluster Formation**

#### "Stable Ph Method"

> Ph Range: 2-12

## **Atomically Precise Metal Cluster:**

> Triphenylphosphine (Tpp) Protected Ag18h16(tpp)2+ Silver Cluster.

### **Proteins:**

> Bovine Serum Albumin (Bsa), Human Serum Albumin (Hsa), Lactoferrin (Lf).

### **Formed Luminescent Metal Cluster:**

> Ag18@bsa And Ag18@hsa.

## **Lactoferrin Stabilizes Smaller Cluster:**

➤ Aq13@lf.

### **Metal Cluster Disintegration**

- Phosphine Protection.
- > Slow Disintegration.

Stable For Months - When Stored At Room Temperature.


### **Enhanced Bioactivity.**

Metal Encapsulated Protein **Exhibits** 

#### **Excellent Luminescent Probe**

For Cancer Cell Imaging.

Confocal microscopic images of internalization of ag18@bsa cluster



# **CONTACT US**

Dr. Dara Ajay, Head TTO Technology Transfer Office, IPM Cell- IC&SR, IIT Madras

# **IITM TTO Website:**

https://ipm.icsr.in/ipm/

Email: headtto-icsr@icsrpis.iitm.ac.in

tto-mktg@icsrpis.iitm.ac.in

Phone: +91-44-2257 9756/ 9719