



### Industrial Consultancy & Sponsored Research (IC&SR)

#### 3D PRINTED THREE-LAYERED POLYMER SCAFFOLD FOR PERIODONTAL REGENERATION, METHOD FOR PREPARING THE SCAFFOLD

IITM Technology Available for Licensing

#### PROBLEM STATEMENT


- **Periodontitis** is a prevalent non-communicable disease affecting teeth, cementum, bone, and periodontal ligament.
- **Current treatments** include plaque elimination, controlling inflammation, and periodontal surgeries.
- **Regeneration of periodontal structure** requires quality regeneration of all three tissues, in the right quantity, and in the same architecture as native tissue.
- **3D printing technology** has gained popularity in **dentistry** for creating **multi-layered scaffolds**, but it does not fully mimic the native periodontal complex and does not guide fiber orientation or anchoring.
- **Various designs** have been proposed for regeneration.

#### Research Lab

Prof. Tuhin Subhra Santra,  
Dept. of Engineering Design

#### TECHNOLOGY

Fig 1 shows a 3D model depicting the different layers in scaffold (a) cross section of the scaffold (b) second scaffold (bone) layer and (c) tri-layer structure



#### TECHNOLOGY CATEGORY MARKET

**Technology:** 3d Printed Three-layered polymer scaffold for periodontal regeneration, method for preparing the scaffold

**Category:** Biotechnology & Genetic Engineering

**Industry:** Biomedical

**Application:** Osteochondral tissue engineering to facilitate multi tissue regenerations.

**Market:** The global market size is **USD 1.55 trillion in 2023** and is projected to **grow at a CAGR of 13.96%** from **2024 to 2030**.

#### INTELLECTUAL PROPERTY

IITM IDF Ref. 2394

Patent No: IN 552910

#### TRL (Technology Readiness Level)

TRL- 3, Experimental proof of concept;

**3D Printed Three-Layer Polymer Scaffold**  
• First scaffold: Vertical polymer scaffold containing angular tubes.

**Second scaffold: Continuous slots to glide over other end of angular tubes.**

**Angular tubes have wide pores for fibroblast migration.**

**Both scaffolds contain polymers from polylactic acid (PLA) and polyethylene terephthalate glycol (PET-G).**

**Printing first and second scaffolds to a platform temperature of 550C to 700C.**

#### CONTACT US


Dr. Dara Ajay, Head TTO  
Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

**IITM TTO Website:**  
<https://ipm.icsr.in/ipm/>

**Email:** [headtto-icsr@icsrpis.iitm.ac.in](mailto:headtto-icsr@icsrpis.iitm.ac.in)  
[tto-mktg@icsrpis.iitm.ac.in](mailto:tto-mktg@icsrpis.iitm.ac.in)

**Phone:** +91-44-2257 9756/ 9719

Figure 2 depicts fluorescence images showing L929 cells growing on PLA scaffold



### Key Features / Value Proposition

#### Layer arrangement:

- First, Second, Third layers.

#### First scaffold:

- Vertical scaffold with angular tubes.
- Vertical scaffold: 6mm x 6mm x 0.3mm to 12mm x 12mm x 0.6mm.
- Polymer weight: First scaffold: 255mg to 265mg;

#### Second scaffold:

- Continuous slots to glide over the first scaffold's angular tubes.
- Polymer weight: Second scaffold: 75mg to 85mg.
- Second scaffold: 6.5mm x 6mm x 0.3mm to 13mm x 12mm x 0.6mm.

#### Angular tubes:

- Wide pores for fibroblast migration and population.
- Pore diameters ranging from 0.1mm to 0.2mm.
- Equal spacing of angular tubes on vertical scaffold.
- Angles: 55°-70°, 90°, 100°-120°
- Angled tubes: 2mm to 4mm, outer diameter 1mm to 2mm, inner diameter 0.8 mm to 1.6mm.

#### Coating:

- Collagen, poly-L-lysine, hydrogel gelatin methacryloyl (GelMA), and fibronectin.

#### 3D Printing Method

- Step (a) printed using a 3D printer at nozzle temperature between 1950°C and 2400°C.
- Step includes coating scaffold with components from collagen, poly-L-lysine, hydrogel gelatin methacryloyl (GelMA), and fibronectin.

#### Technique used

- Fused deposition modeling (FDM).

#### Coating Scaffold with Components

- Collagen, poly-L-lysine, Hydrogel Gelatin Methacryloyl (GelMA), and fibronectin.

### CONTACT US

Dr. Dara Ajay, Head TTO  
Technology Transfer Office,  
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:  
<https://ipm.icsr.in/ipm/>

Email: [headtto-icsr@icsrpis.iitm.ac.in](mailto:headtto-icsr@icsrpis.iitm.ac.in)  
[tto-mktg@icsrpis.iitm.ac.in](mailto:tto-mktg@icsrpis.iitm.ac.in)

Phone: +91-44-2257 9756/ 9719