

IIT MADRAS Technology Transfer Office TTO - IPM Cell



Industrial Consultancy & Sponsored Research (IC&SR)

A SCAFFOLD OF NANOFIBROUS MATRIX AS A DIABETIC WOUND DRESSING

IITM Technology Available for Licensing

Problem Statement

Indian Institute of Technology Madras

- In the present subject matter, the **problem** in electrospinning agarosecurdlan blend is that agarose & curdlan are soluble in hot aqueous solution at higher temperature, however they gel/solidify at room temperature and hence such a blend could not be employed for the production of nanofibers by means of electrospinning.
- The present invention addresses said problem and provides a means for the production of scaffold.

Technology Category/Market

Technology: Scaffold of nanofibrous matrix; Industry & Application Biotechnology & pharmaceutical companies; Biomedical Device; Market: The global scaffold technology market is projected to grow USD 5.09Bn by 2031, at a CAGR of 18% during the forecast period (2024-31)

Technology

- Present patent is related to a scaffold comprising a matrix of nanofibers, which is meant for healing chronic wounds and is prepared from a blend of biocompatible polymers using an **electrospinning** technique.
- Further, said patent explains a method for producing a scaffold for chronic wound healing.
- Present Patent comprises a **matrix** of nanofibers of biocompatible polymers comprising agarose, curdlan, polyvinyl alcohol with or without antibacterial drug.
- Said **nanofibers** are a blend of biocompatible polymers comprising
- 1%-10% w/v of agarose in formic acid or DMSO;
- 1%-5% w/v of curdlan in formic acid or DMSO; and

CONTACT US

Dr. Dara Ajay, Head TTO Technology Transfer Office, IPM Cell- IC&SR, IIT Madras

IITM TTO Website: https://ipm.icsr.in/ipm/

- 10%-12% w/v of polyvinyl alcohol in water,
- wherein the ratio of combination of agarose in formic acid or DMSO and curdlan in formic acid or DMSO to polyvinyl alcohol in water in the blend ranges from 1:2 to 1:4.



Fig.1 illustrates the scanning electron microscopy images (Scale- 1µm) (a) ACP, ACP1%CFX, (c) ACP3%CFX, (d) (b) ACP5%CFX

Intellectual Property

IITM IDF Ref. 2427; IN Patent No. 523793 (Granted)

TRL (Technology Readiness Level)

TRL-3, Experimental proof of concept;

Research Lab

Prof. Sampath Kumar T S; Dept. of Metallurgical and Materials Engg. **Prof. Mukesh Doble;** Dept. of Biotechnology

> Email: smipm-icsr@icsrpis.iitm.ac.in sm-marketing@imail.iitm.ac.in Phone: +91-44-2257 9756/ 9719



Indian Institute of Technology Madras



Industrial Consultancy & Sponsored Research (IC&SR)





FIG.2: Depicts the Disc diffusion assay showing antibacterial efficacy of fabricated scaffold against E. coli

(a, b, c, and d for ACP, ACP1%CFX, ACP3%CFX, & ACP5%CFX respectively) and S. aureus(e, f, g, and h for ACP, ACP1%CFX, ACP3%CFX, and ACP5%CFX respectively).

FIG.3a & 3b: Depicts In vitro gene expression levels of

- (a) Inflammatory cytokines,
- (b)Anti inflammatory cytokines, in THP-1 macrophages on treatment with ACP, and ACP3%

Key Features / Value Proposition

* <u>Technical Perspective</u>:

Important Features:

- The nanofibers of biocompatible polymers further comprise antibacterial drug, and it helps combat infection and hasten the healing cascade.
- Said scaffold effectively absorbs the exude emanating from chronic wounds such as fungating wounds, fistulae, venous leg ulcers, and diabetic foot ulcers, and poses no exude leakage problem, & specifically advantageous for healing said chronic wounds.
- The scaffold facilitates easy transmission into spacerestricted wounds (cavity wounds) with areater flexibility.
- The matrix of **nonwoven nanofibers** have diameter in the range of 115-146 nm with effective exude absorption ability (~ 450 - 500 % swelling).
- Exemplary: Cell viability was ~90% using mouse fibroblasts (L929 cell line) i.e. better cell viability (Refer Fig 4a)



Fig. 4a: Shown Percent metabolic activity of fabricated scaffold with (a) L929 mouse fibroblasts:

CONTACT US

Dr. Dara Ajay, Head TTO Technology Transfer Office, IPM Cell- IC&SR, IIT Madras

IITM TTO Website: https://ipm.icsr.in/ipm/ Email: smipm-icsr@icsrpis.iitm.ac.in sm-marketing@imail.iitm.ac.in Phone: +91-44-2257 9756/ 9719