



GLYMO FUNCTIONALIZED SILICON NANOPOROUS MEMBRANE (SNM) FOR NON-INFLAMMATORY, NON-IMMUNOGENIC HEMODIALYSIS APPLICATIONS AND METHOD THEREOF

IITM Technology Available for Licensing

Problem Statement

- **Hemodialysis membranes** help in uremic toxin clearance **for kidney treatment**.
- **Conventional polymer-based** membranes are thick, inconsistent, and **prone to fouling**, leading to **inefficient toxin clearance** and potential **blood protein loss**.
- Further, traditional **silicon-based membranes** with PEG coating suffer from **high costs**, **instability** under blood exposure, and **inadequate precise pore size control**.
- There is a **need for thin, stable and reusable membrane** that **minimizes fouling** while ensuring **toxin clearance** and **protein retention**.

Intellectual Property

- IITM IDF Ref 2553
- IN 559892 Patent Granted

TRL (Technology Readiness Level)

TRL 4 Technology validated in Lab

Technology Category/ Market

Category- Drugs and Pharmaceutical Engineering

Industry Classification:

Healthcare and Medical Devices Industry;
Nanotechnology and Advanced Materials;
Advanced membrane and filtration technologies

Applications:

Hemodialysis devices; Blood purification systems; Modular and scalable dialysis apparatuses

Market report:

The global hemodialysis membrane market was valued at USD 9.16 billion in 2024 and is projected to grow to USD 12.93 billion by 2031, with a CAGR of 5.1%

Research Lab

Prof. Enakshi Bhattacharya
Dept. of Electrical Engineering




Figure: Surface functionalization of SNM with GLYMO. This procedure completely reduced the surface interaction with major uremic toxins, particularly urea, and proteins.




Figure: Optical image of single SNM (left) and magnified TEM image of SNM (right). The small pore-size and high porosity helps in efficient clearance of small and medium sized biomolecules (MW < 65 kDa) at a faster rate, while retaining the higher MW proteins in blood




Figure: SNM array: 6x6 array (left) and 5x5 array (right). This greatly improved the exchange area and surface area to volume ratio for dialysis, thereby reducing the clearance time for higher volumes

CONTACT US

Dr. Dara Ajay, Head TTO
Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:
<https://ipm.icsr.in/ipm/>

Email: headtto-icsr@icsrpis.iitm.ac.in
ttooffice@icsrpis.iitm.ac.in

Phone: +91-44-2257 9756/ 9845




Figure: ATR comparison of GLYMO treated membrane surface before and after the exposure to the human blood. The stability was checked by the presence of Si-O-Si peak at $\sim 1200 \text{ cm}^{-1}$ in FTIR, which is the characteristic indication of silane ether bond formed after functionalization. Histological blood slides were prepared with the blood sample after the test and screened under a microscope. The cells were intact indicating that there was no chemical effect or shear stress of the functionalized membranes on the blood cells.



Figure: Irritability, pyrogenicity and systemic toxicity of the GLYMO treated Silicon membranes were tested in Lister mice. As per standard procedure, the animals were kept for 7 days in quarantine before testing. In test animal, a cut was made on the skin, empty space was created under the skin and 2 cm^2 area of functionalized silicon membranes was inserted. The histological analysis of the surrounding tissue and the organ tissue showed no inflammatory or sensitization response

Technology

The invention is a GLYMO-functionalized Silicon Nanoporous Membrane (SNM) designed for hemodialysis, ensuring efficient uremic toxin clearance while retaining vital blood proteins.

Fabricated via ICP-CVD and rapid thermal processing, the SNM is 15 nm thick with $\sim 10 \text{ nm}$ pore diameters, offering high precision and scalability.

GLYMO silanization converts surface groups to stable Si-O-Si bonds, minimizing biofouling and enhancing long-term stability under blood exposure.

The membrane clears 50% of uremic toxins in 15 minutes, meets ISO 10993-4:2017(E) standards, and demonstrates excellent hemocompatibility with minimal cytotoxicity.

Targeted for advanced hemodialysis devices, the technology offers reusability, durability, and a compact design, providing a safer and more efficient solution for kidney failure treatment.

Key Features / Value Proposition

- The invention uses a 15 nm-thin SNM with $\sim 10 \text{ nm}$ pores, ensuring faster uremic toxin clearance and reduced clogging compared to conventional polymer membranes that are thicker and less uniform.
- GLYMO functionalization minimizes biofouling by stabilizing the surface with Si-O-Si bonds, enhancing biocompatibility and reducing adverse blood interactions.
- Compare to conventional polymer or PEG-coated membranes, the invention efficiently retains vital blood proteins while clearing toxins, lowering risks like sepsis and hemolysis.
- The fabrication process in the invention using ICP-CVD and rapid thermal processing is scalable and more cost-effective than expensive, complex nanolithography techniques.
- Rigorous in vitro tests confirm the invented membrane's non-inflammatory, non-immunogenic performance, demonstrating superior durability and reusability for safe, long-term hemodialysis applications.

CONTACT US

Dr. Dara Ajay, Head TTO
Technology Transfer Office,
IPM Cell- IC&SR, IIT Madras

IITM TTO Website:
<https://ipm.icsr.in/ipm/>

Email: headtto-icsr@icsrpiis.iitm.ac.in
ttooffice@icsrpiis.iitm.ac.in

Phone: +91-44-2257 9756/ 9845